

# Medik8<sup>®</sup>

SKINCARE BIOTECHNOLOGY™

**Titanium Gold  
Microneedle  
Roller System**

Restore a healthy-looking complexion with painless microneedle therapy. Stimulate natural collagen formation to reduce the appearance of lines and wrinkles. Encourage cell repair pathways to improve the appearance of scars, especially indented acne scars. Improve penetration and performance of Medik8 serums and creams as part of your Medik8 skincare regime.

### Simple Directions For Use

(for full directions go to [www.medik8.com](http://www.medik8.com))

1. Sterilise the Medik8 roller with a sterilising solution and allow to dry before every use.
2. Cleanse face, neck, décolleté or body using a Medik8 cleanser suitable for your skin type and pat skin dry.
3. Using a constant, gentle rolling motion roll 10 times over each skin area.
4. Apply your preferred Medik8 Prevent and/or Correct product, following the product instructions.
5. Apply Hydr8 SPF25 or a suitable sunscreen.
6. Clean roller with hot water and return to case, dry.

At first repeat treatment twice a week. Before increasing frequency of use, consult a professional therapist .

**Replace roller after 3 months of use.**

### Cautions

Following treatment, the skin will appear flushed. This should subside within a few hours. Cease use if you develop a skin rash or irritation and resume treatment after skin feels comfortable. Do not use on the following skin conditions: active acne, eczema, rosacea, psoriasis, actinic (solar) keratosis and raised moles or warts. Do not use if pregnant. Do not use if on blood thinning medication or if you are a haemophiliac. Store out of reach of children. Not to be used by anyone under the age of 18.

### Manufactured by

Pangaea Laboratories Ltd.,  
58 Crewys Road,  
London NW2 2AD  
Tel: +44 (0)845 673 2222  
[www.medik8.com](http://www.medik8.com)

# nanogen<sup>®</sup>

## Microneedle Scalp Roller System

Greatly increase the absorption of topical serums like nanogaine and minoxidil and through massage promote increased blood flow to the scalp for healthy hair growth. There is also some evidence that causing minor 'trauma' to scalp skin through micro-needling can also promote an automatic repair response which is positive for hair follicle growth.

### Directions for Use

1. For best results, cleanse the scalp removing oil, dirt, sebum and flakes using Onescalp. Rinse thoroughly and towel dry.
2. Use a constant, gentle rolling action. Move in one direction only, from the crown to front and /or sides.
3. Apply topical solutions like nanogaine or minoxidil, gently massaging into the scalp following the product instructions. Use roller at night as generally topical products require 4 hours to absorb fully.
4. Ensure roller is cleaned thoroughly after each use. Run under hot or near boiling water for hygienic purposes and allow to dry before storing in case. Always allow to cool to room temperature before use on skin.

### Cautions

Do not share your Nanogen Scalp Roller with anyone else.

Do not use more than once daily.

As a precaution, do not use if pregnant or if taking any blood thinning medication.

Store in a safe place out of reach of children.

Not to be used by persons under the age of 18.

### Manufactured by

Pangaea Laboratories Ltd,  
58 Crewys Road,  
London NW2 2AD  
Tel: +44 (0)845 673 2222  
[www.nanogen.co.uk](http://www.nanogen.co.uk)